



Mixing



Labfirst Scientific Instruments (Shanghai) Co., Ltd.

Vortex C

Cell Disruptor

Specifications	Vortex C
Functions	Disrupts cells, yeast, bacteria, algae, fungi, etc. to prepare cell lysates for further processing
Mixing motion	High-speed circumferential vortex and oscillation function
Orbital diameter	4mm
Speed range	0-3000rpm (Adjustable)
Speed display	Scale
Capacity	8×2 mL
Power Requirements	200-240V/ 100-120V, 50/60 Hz, 60W
Dimension	127X130X160mm
Weight	2.0kg



Features

- Durable, reliable and affordable
- Easy to operate, rapid hands-free disruption of up to 8 microtubes
- Processes at a constant, controlled shear rate
- Maintains genome integrity of the biological samples
- Facilitates subsequent purification of genomic material
- Best suited for research laboratories demanding accurate, reproducible, repeatable results

Applications

Suitable for all sample types; Lysis of yeast and fungi; for Isolation of genomic material from faecal samples; fungi and bacteria samples; soil samples; for Isolation of PCR-Ready Genomic DNA from Soil Samples; DNA & RNA Purification; Lysis of bacteria; Lysis of soft tissues; Isolation of genomic DNA; Purification of DNA; Cell Lysis; Used to lyse cell wall of yeast; Break hard/brittle samples.

VT1.3.1 tube adapter

Cat. No.921004

for 48 test tubes, Ø6mm (1 ml)


VT1.2 Tube holding rod

Cat. No.921002

used with tube adapters


VT1.3.2 tube adapter

Cat. No.921005

for 15 test tubes, Ø10mm (1.5-2 mL)


VT1.3 Universal top plate

Cat. No.921003

Ø 100mm


VT1.3.3 tube adapter

Cat. No.921006

for 16 test tubes, Ø12mm (5 mL)


VT1.3.6 Platform pad

Cat. No.921009

Ø99mm, multiple small tubes and vessels


VT1.3.4 tube adapter

Cat. No.921007

for 8 test tubes, Ø16mm (10 mL)


VT1.3.7 vacuum chuck

Cat. No.921010

made of rubber


VT1.3.5 tube adapter

Cat. No.921008

for 8 test tubes, Ø20mm (10-15 mL)


PS1.2 Microplate clamp

Cat. No.921202


VT1.3.6 PCR tube adapter

Cat. No.921012

 for 56 test PCR tubes, strips of 8 ,
Ø5mm (0.1mL-0.5 mL)


Application of accessories

Accessories	Adjustable speed model (0-2500rpm)	
	Touch mode (High speed area)	Continuous mode (Low speed area)
VTI.1	Y	Y
VTI.2 + VT1.3.1	Y	
VTI.2 + VT1.3.2	Y	
VTI.2 + VT1.3.3	Y	
VTI.2 + VT1.3.4	Y	
VTI.2 + VT1.3.5	Y	
VTI.3 + VT1.3.6	Y	
VTI.3 + VT1.3.1		Y
VTI.3 + VT1.3.2		Y
VTI.3 + VT1.3.3		Y
VTI.3 + VT1.3.4		Y
VTI.3 + VT1.3.5		Y



Lab1st mixers range can be divided into different series as per mixing speed, working mode, load capacity etc:

Series	Image	Model	Speed Range	Working mode	Capacity
Vortex Mixer		Vortex S	0~3000rpm	Orbital (touch + continuous modes)	-
		Vortex F	3000rpm(fixed)	Orbital (touch + continuous modes)	-
		Vortex E	3000rpm (fixed)	Orbital (touch mode)	-
		Vortex C	0~3000rpm	Orbital (touch + continuous modes)	8x2mL
Microplate Mixer		Microster M	0~1500rpm	Orbital (continuous mode)	0.5Kg
Disk Rotator		Rotaster 3	10~70rpm	360° rolling (timer + continuous)	60x1.5mL, 16x15mL, 8x50mL, (available for 3 layers)
Rotisserie Rotator		Grillster T	10~70rpm	360° rolling (timer + continuous)	48x1.5mL, 24x15mL, 24x50mL
Disk Rotator		Rotaster 1 basic	0~80rpm	360° rolling (continuous)	48x2mL, 16x15mL, 8x50mL
Rotisserie Rotator		Grillster basic	0~80rpm	360° rolling (continuous)	32x2mL, 20x15mL, 16x50mL
Rocker		See-saw 30 basic	0~30rpm	30±3° rocking (continuous)	10x10mL, 4x50mL
		See-saw 30L	0~30rpm	30±3° rocking (continuous)	16x10mL 6x50mL
		See-saw 30D	0~30rpm	40±3° rocking (continuous)	32x10mL, 12x50mL
Roller		Rollster 6	0~70rpm	360° rocking&rolling (timer + continuous)	4Kg
		Rollster 6 basic	0~70rpm	360° rocking&rolling (continuous)	4Kg

Lab1st shaker series can be divided to different types as per different motion models:

Series	Image	Model	Plate Size	Display	Speed Range	Max. Load Capacity
Orbital Shaker		Orbister 330	32x32cm	LCD	100-500rpm	7.5Kg
		Orbister 250	25x24.5cm	LCD	100-800rpm	2.5Kg
		Orbister 330 basic	33.5x33.5cm	LCD	70-400rpm	3Kg
		Orbister 250 basic	26.8X26.8cm	LED	40-200rpm	2Kg
Linear Shaker		Recipster 330	32x32cm	LCD	100-350rpm	7.5Kg
		Recipster 250	25x24.5cm	LCD	100-350rpm	2.5Kg
Rocking Shaker		Wavester 330	30x28.5cm	LCD	10-70rpm	10Kg
		Wavester 250 basic	26.8x26.8cm	LED	10-80rpm	2Kg
3D Shaker		3D-Ster 330	30x28.5cm	LCD	10-70rpm	5Kg